## IN THE CLAIMS

Please amend the claims as follows:

Claims 1-18 (Canceled).

Claim 19 (Currently Amended): A <u>blood analyzer for analyzing whole blood which</u> includes a quality control device for a blood analyzer using whole blood, <u>said blood analyzer</u> comprising:

means for analyzing patient blood samples;

a quality control device including,

means for storing by refrigeration [[for]] control <u>bloods</u>, <del>bloods</del>;

means for restoring to temperature the control bloods to a temperature prescribed by [[the]] <u>a</u> manufacturer of the control <u>bloods</u>, <del>bloods</del>;

means for stirring <u>the control bloods</u> for re-suspension of <u>cells</u>, <del>cells</del>; and means for sampling the <u>control bloods</u> <del>blood thus prepared</del>,

wherein said quality control device is incorporated into the blood analyzer.

Claim 20 (Currently Amended): A <u>blood analyzer device</u> according to claim 19, wherein the means for storing control bloods comprises a specified number of tubes sealed by a bung and arranged in a tube support in contact with a refrigeration block configured to adjust temperature and maintain an optimum temperature for storing the control bloods.

Claim 21 (Currently Amended): A <u>blood analyzer</u> device according to claim 20, wherein the refrigeration block is a Peltier effect refrigeration block.

Claim 22 (Currently Amended): A <u>blood analyzer device</u> according to claim 20, wherein the <u>tube support is disconnected from the refrigeration block for means for restoring</u> the temperature of the control bloods <u>includes the tube support configured to disconnect from the refrigeration block</u>.

Claim 23 (Currently Amended): A <u>blood analyzer device</u> according to claim 21, wherein the means for restoring includes the Peltier effect refrigeration block configured to receive current which supplying the Peltier effect refrigeration block is interrupted for a specified period of time such that for restoration of the temperature of the control bloods is restored.

Claim 24 (Currently Amended): A <u>blood analyzer</u> device according to claim 21, wherein the Peltier effect refrigeration block is controlled to <u>be</u> reset and maintain quality control to [[its]] a utilization temperature according to specifications of the manufacturer.

Claim 25 (Currently Amended): A <u>blood analyzer device</u> according to claim 20, wherein the means for [[the]] stirring <u>includes the tube support articulated about a hinge of the refrigeration block and configured to operate by inverting operates by rocking and/or inversion formed by the tube support articulated about a hinge of the refrigeration block.</u>

Claim 26 (Currently Amended): A <u>blood analyzer device</u> according to claim 25, wherein an angle of inversion of the <u>tube support</u> is between 100° and 180°.

Claim 27 (Currently Amended): A <u>blood analyzer</u> device according to claim 19, wherein the means for <u>storing stirring</u> includes low-speed Vortex stirring means.

Claim 28 (Currently Amended): A <u>blood analyzer</u> device according to claim 19, wherein the means for sampling includes a needle configured to draw blood from [[the]] tubes.

Claim 29 (Currently Amended): A <u>blood analyzer device</u> according to claim 28, wherein the needle is <u>configured to be</u> driven in a transverse movement over [[the]] tubes of <u>patient</u> blood <u>samples</u> to be analyzed and the control bloods as well as over a counting block comprising mixing and rinsing tanks and <u>is configured to be driven</u> in a vertical movement to penetrate into the tubes by piercing [[the]] bungs or by descending into the counting block comprising mixing and rinsing tanks to carry out rinsing or dilutions of the blood.

Claim 30 (Currently Amended): A <u>blood analyzer</u> device according to claim 29, wherein <u>the bungs are configured to be pierced</u> the <u>piercing of the bungs is effected</u> when the tubes on their a support are in a high or low position.

Claim 31 (Currently Amended): A <u>blood analyzer device</u> according to claim 19, further comprising programmable processing means for checking that values obtained by passing through <u>each a quality control procedure</u> correspond to limit values and expected values of the control blood.

Claim 32 (Currently Amended): A <u>blood analyzer</u> device according to claim 31, wherein the processing means triggers an alarm [[if]] <u>when</u> the values obtained during running of the quality control <u>procedure</u> are outside the expected <u>values</u> <u>limits</u>.

Claim 33 (Currently Amended): A <u>blood analyzer device</u> according to claim <u>37</u> [[19]], <u>further comprising means for wherein the</u> triggering <u>unit a quality control procedure</u> triggers the quality control procedure either directly by an operator <u>input</u>, [[or]] automatically, or via an external connection to a control unit.

Claim 34 (Currently Amended): A <u>blood analyzer</u> device according to claim 19, wherein transfer and analysis of data are affected via an internal or external network implementing standards of HL7, ASTM, or XML.

Claim 35 (Currently Amended): A <u>blood analyzer device</u> according to claim 20, wherein the tubes include means for identification and tracking by barcodes, electronic chips, and/or or magnetic labels for identifying and tracking the tubes.

Claim 36 (Canceled).

Claim 37 (New): A blood analyzer according to claim 19, wherein the quality control device further includes a triggering unit which triggers a quality control procedure to determine whether the analyzer is functioning properly based on a comparison using the control bloods.

Claim 38 (New): A blood analyzer for analyzing whole blood which includes a quality control device, said blood analyzer comprising:

an analyzing unit which analyzes patient blood samples;

a quality control device including,

a storing unit which stores, by refrigeration, control bloods,

Application No. 10/564,784 Reply to Office Action of May 11, 2010

a restoring unit which restores the control bloods to a temperature prescribed by a manufacturer of the control bloods,

a stirring unit which stirs the control bloods for re-suspension of cells, and a sampling unit which samples the control bloods.